

PMT[®] HALO SYSTEMS

C E R V I C A L T R A C T I O N



Model 1233
w/ Lambs Wool Liner



Model 1233
w/ Coolmax[®] Liner

The PMT[®] Corporation's Halo system is the most reliable, easy to use and effective Halo structure available. This Halo system is designed to provide optimal patient comfort and healing. Compared to spinal fusion, Halo fixation allows the body to naturally repair itself at a reduced cost.

The PMT Halo head ring utilizes hand wound graphite to ensure product quality and consistency. The proprietary titanium skull pins are specially designed with a parabolic tip to increase the anchoring qualities and reduces the need to retighten the pins as with other designs. The medical grade titanium, TI 6AL-4V ELI, combined with the graphite head ring and support posts are MRI conditional at both 1.5T and 3.0T as well as CT compatible to provide the best imaging characteristics.

PMT[®] CORPORATION
QUALITY, INNOVATION AND CUSTOMER DRIVEN